



HEALTHCARE COMMUNICATION SOLUTIONS, INC.

610-743-4668 • info@hcsinc-pa.com
www.hcsinc-pa.com

2669 Shillington Road • Suite 331 • Sinking Springs, PA 19608



The CTD Shield is a single use disposable product that can help reduce cross-contamination risks associated with in-room hand held devices. The CTD Shield effectively reduces or eliminates repairs to the patient control units due to liquids, bodily fluids and chemicals when properly installed. In addition, the CTD Shield will help reduce cleaning and replacement expense, as well as, provide peace of mind for the patients and staff. The CTD Shield is constructed from 1.5 mil polyurethane film. The product is approximately 4 1/2 inches wide and 12 inches long.

POLYURETHANE PROTECTIVE BARRIER Designed for Hand Held Nurse Call Devices

The CTD Shield is designed to be used with Gen II & Gen III Curbell Pillow Speakers. The primary use is to help facilitate contact transmitted infection control by reducing the risks associated with bacterial, viral, chemical, and other infectious and hazardous contaminants. The CTD Shield also provides a fluid barrier to prevent fluids from entering the porous surface of interior of the control unit which should help extend the life of the unit. For optimum protection against contact transmitted diseases, the shield should be changed after every patient discharge.

